

SolaCare Psychiatric Symmetry



MEDTRICA

Multi Layered Foam Support

Top foam

1.9 lb density and IFD 27-35. Mattress is Reversible.

Perimeter Foam Support

2.0 lb density and IFD 86-94.

AVAILABLE SIZES:


84" x 35" x 6"

80" x 35" x 6"

78" x 35" x 6"

76" x 35" x 6"

Custom Sizes



Fluid Proof Cover Options:

- Dartex™ 4 way stretch
- 70D Nylon

Base Foam Support

1.8 lb density and IFD 36-44. Provides a firm support base that prevents "bottoming out" even as the head section is fully elevated.

NON-ALLERGENIC

The entire foam core is constructed of 100% medical grade Green Friendly non-allergenic foam.

ANTI-MICROBIAL

All components are anti-microbial and assist in prohibiting growth of bacteria and fungi.

FLAMMABILITY

Meets Federal Standard 16 CFR 1632, Cal. TB 117. Optional barrier to meet Cal. TB 129 and Boston Fire IX-11.

COVER DESIGN

Fully enclosed cover, No Handles, No Flaps, No Zippers.

STANDARD MATTRESS COVER

Durable, snug fitting low denier 70D Nylon is soft, moisture vapor permeable, and resistant to staining with low slip vinyl bottom cover. Will not allow fluids to pass through and can be easily cleaned with standard hospital cleaning products. 100% 70D Nylon, no butyl backing.

OPTIONAL MATTRESS COVER

Dartex™ 4 way Stretch Pressure Reducing fabrics are the world leader in healthcare fabrics. Dartex™ fabrics are engineered to provide the best in MVT, pressure management, abrasion and tear resistance while maintaining optimal performance throughout the entire lifespan of the product.

GUARANTEE

A full non prorated guarantee on all SolaCare Series mattresses instead of the industry standard limited warranty.

- 3 year cover guarantee
- 3 year foam core guarantee

SAFETY, COMFORT & DURABILITY

The SolaCare Psychiatric Symmetry Mattress is specifically engineered to provide lasting optimum performance w/ the highest quality pressure reducing materials.

Comfort, Safety, Durability and Budget is addressed. Fluid-proof Dartex™ or 70D Nylon fabric covers.

info@medtrica.com
www.medtrica.com



A better quality mattress, more durable, safe and comfortable - for a better price.