


Standard Stretcher Mattress



MEDTRICA

Fluid Proof Cover Options:

- Dartex™ 4 way stretch
- Vinyl



Anti-Shear Liner Option:

Dissipates Heat and Reduces Friction

Foam Core

- 3", 4" or 5" thick
- 1540 foam

AVAILABLE SIZES:

All standard and Custom sizes Available to Fit all Makes and Models of Stretcher Frames

GUARANTEE:

A full non-prorated guarantee on both the cover and the core for all stretcher mattresses.

Adhesive Velcro included to keep the mattress secure to the stretcher.

FLAMMABILITY

All materials meet Cal TB117 and CFR1632. Option to meet Cal TB129 & Boston FC IX 11.

All Stretcher mattresses meet and exceed OEM specifications.

CHOICE OF MATTRESS CONSTRUCTION:

- 3", 4" or 5" Thick; all standard and custom Lengths & Widths
- Tapered or Square corners.

COVERS:

All covers are 100% medical grade fabric-anti-fungal, anti-bacterial, latex free and completely fluid-proof. Choices are a durable Vinyl or a Pressure Reducing stretchable Dartex™ fabric.

DARTEX™ 4 WAY STRETCHABLE PRESSURE REDUCING FABRICS:

Dartex™ 4 way Stretchable Pressure Reducing Fabrics are the world leader in Healthcare cover fabrics. Dartex™ fabrics are engineered to provide the best in MVT, pressure management, abrasion and tear resistance while maintaining optimal performance throughout the entire lifespan of the stretcher mattress.

FOAM CORES:

100% medical grade Green Friendly non-recycled, non-allergenic, HD polyurethane foam without fillers or additives.

CLEANING:

Stretcher mattresses are disinfected and wiped clean using standard cleaning products.

MORE DURABLE & VERSATILE

Medtrica stretcher mattresses provide greater comfort & durability because of unsurpassed Quality.

100% medical grade Green Friendly foam core.

Manufactured in the USA & Canada - Engineered to Last.

info@medtrica.com
www.medtrica.com



A better quality stretcher mattress, more durable - for a better price.