





APPLI-KIT™ BEDSIDE PRECLEAN KIT

Description	Code	
<p>For optimal rinsing and pre-cleaning of instruments at bedside</p> <p>Delivers Low Foaming D-Zyme™ enzymatic detergents containing 3 enzymes - Protease, Amylase and Lipase, PLUS Surfactant Detergents to break down proteinaceous bioburden during the preclean of surgical instruments</p> <ul style="list-style-type: none"> • Green Friendly, 100% Medical Grade foam sponges, No Flaking or Scratching • Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacteriostatic <p>Stable wide base pouch for easy fluid transfer. Clear pouch bag for a direct view of the CONTENTS</p> <p>Add up to 500 ml of water to the POUCH then soak and squeeze the sponge 3 - 5 times to activate the D-Zyme™ enzymatic concentrate.</p>		 <p>Appli-Kit™ Flat- Add Water to bag</p>
<p>Appli-Kit™ w/ D-Zyme™ Concentrate in stand-up pouch bag</p> <p>Options:</p>		
<ul style="list-style-type: none"> • Appli-Kit™ Flat w/ D-Zyme™ Concentrate - Add Water to pouch (up to 500 ml) 100 ct / 1 case 	<p>AFC500100</p>	
<ul style="list-style-type: none"> • Appli-Kit™ Tubular w/ D-Zyme™ Concentrate - Add Water to pouch (up to 500 ml) 70 ct / 1 case 	<p>ATC500070</p>	 <p>Appli-Kit™ Tube - Add Water to bag</p>

Manufactured by:

MEDTRICA SOLUTIONS LTD.
104A-6350 204th Street
Langley, BC Canada V2Y 2V1
t: 604-618-2291



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
101-128
Redondo Beach, CA 90278
t: 877-703-1079