



Medtrica Overlays - O.R.



	Convuluted Foam Overlay	Flat Foam Overlay	Convuluted Memory Foam	Flat Memory Foam
				
	Convuluted memory foam - uncovered	Flat Foam Overlay - uncovered	Convuluted memory foam - uncovered	Flat Memory Foam - covered
Cover Options	<ul style="list-style-type: none"> •Uncovered •70D Nylon •Dartex 4-way stretch 	<ul style="list-style-type: none"> •Uncovered •70D Nylon •Dartex 4-way stretch 	<ul style="list-style-type: none"> •Uncovered •70D Nylon •Dartex 4-way stretch 	<ul style="list-style-type: none"> •Uncovered •70D Nylon •Dartex 4-way stretch
Foam	High Density FR Rated	High Density FR Rated	2.4 lb or 5 lb Viscoelastic 'Memory' Foam	2.4 lb or 5 lb Viscoelastic 'Memory' Foam
Sizes	8"x8"x2" - pack of 32 10"x8"x0.75" - pack of 48 12"x12"x6" - pack of 8 16"x7.5"x2" - pack of 16 20"x7.5"x2" - pack of 16 20"x12"x2" - pack of 8 21"x18"x7.5" - pack of 2 35"x18"x2.25" - pack of 4 35"x18"x2.25" - pack of 8 72"x20"x2" - pack of 2 20"x10"x6" HD Foam - pack of 4 20"x6"x6" HD Foam - pack of 8	8"x8"x2" - pack of 32 10"x8"x0.75" - pack of 48 12"x12"x6" - pack of 8 16"x7.5"x2" - pack of 16 20"x7.5"x2" - pack of 16 20"x12"x2" - pack of 8 21"x18"x7.5" - pack of 2 35"x18"x2.25" - pack of 4 35"x18"x2.25" - pack of 8 72"x20"x2" - pack of 2 26"x6"x1" - pack of 32 26"x6"x2" - pack of 16 26"x6"x3" - pack of 10 26"x6"x4" - pack of 8	16"x5.5"x2" - pack of 10 20"x5.75"x2.25" - pack of 8 24"x5.75"x2.25" - pack of 8 26"x6"x2" - pack of 8 26"x6"x3" - pack of 8 26"x6"x4" - pack of 4	15"x6.5"x0.25" - pack of 48 16"x5.5"x2" - pack of 10 20"x5.75"x2.25" - pack of 8 24"x5.75"x2.25" - pack of 8 24"x4.5"x0.375" - pack of 32 24"x6.5"x0.375" - pack of 32 26"x6"x1" - pack of 16 26"x6"x2" - pack of 8 26"x6"x3" - pack of 8 26"x6"x4" - pack of 4 46"x20"x0.5" - pack of 4 72"x20"x0.5" - pack of 4
Anti-Bacterial, anti-Fungal, anti-static Stain & Tear resistant, Latex Free	Yes	Yes	Yes	Yes
100% Medical Grade Foam and Fabric	Yes	Yes	Yes	Yes
Anti-Shear Inner Liner	Optional	Optional	Optional	Optional
Low Slip Bottom	Optional	Optional	Optional	Optional
Custom Sizing, Shapes and covers available	Yes	Yes	Yes	Yes
Fire Rating	All materials exceed fire standards Cal TB 117, 16CFR-1632	All materials exceed fire standards Cal TB 117, 16CFR-1632	All materials exceed fire standards Cal TB 117, 16CFR-1632	All materials exceed fire standards Cal TB 117, 16CFR-1632
Additional Fire Barrier	Optional	Optional	Optional	Optional
Warranty	70D Nylon cover - 1 year Dartex cover - 1 year	70D Nylon cover - 1 year Dartex cover - 1 year	2.4lb foam w/70D Nylon or Dartex cover - 1 year 5.0lb foam w/70D Nylon or Dartex cover - 2 year	2.4lb foam w/70D Nylon or Dartex cover - 1 year 5.0lb foam w/70D Nylon or Dartex cover - 2 year