

# Bedside Cleaning Ready to Use



MEDTRICA

## APPLI-KIT™ BEDSIDE PRECLEAN KIT W/ PRE-MIXED D-ZYME SOLUTION

Description	Code	
<p>For optimal rinsing and pre-cleaning of instruments at bedside</p> <p>Delivers Low Foaming d-zyme enzymatic detergents containing <b>3 enzymes</b> - Protease, Amylase and Lipase, <b>PLUS Surfactant Detergents</b> to break down proteinaceous bioburden during the preclean of surgical instruments</p> <ul style="list-style-type: none"> <li>• Green Friendly, 100% Medical Grade foam, No Flaking or Scratching</li> <li>• Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacteriostatic</li> <li>• Pre-mixed d-zyme enzyme detergents solution provided in a stand-up pouch bag</li> <li>• Stable wide clear pouch for easy fluid transfer and direct view of the solution.</li> </ul> <p>The water that is mixed with d-zyme enzymatic detergents concentrate to make the pre-mixed solution goes through a multi-stage cleaning process in the following order: 1) Particulate Filtration; 2) Charcoal Filtration; 3) Semi Permeable Membrane (reverse osmosis); 4) UV Light</p>		 <p>Appli-Kit™ Flat - 500ml pre-mixed</p>
<p><b>Options:</b></p> <p><b>Appli-Kit™ Flat:</b></p> <ul style="list-style-type: none"> <li>• <b>215 ml</b> of pre-mixed d-zyme Enzymatic Detergent Solution - 50 / case</li> <li>• <b>500 ml</b> - 40 / case</li> </ul>		 <p>Appli-Kit™ Tubular - 215ml pre-mixed</p>
<p><b>Appli-Kit™ Tubular:</b></p> <ul style="list-style-type: none"> <li>• <b>215 ml</b> of pre-mixed d-zyme Enzymatic Detergent Solution - 35 / case</li> <li>• <b>500 ml</b> - 25 / case</li> </ul>	<p>APF215050 APF500040</p> <p>APT215035 APT500025</p>	

Manufactured by:

MEDTRICA SOLUTIONS LTD.  
104A-6350 204th Street  
Langley, BC Canada V2Y 2V1  
t: 604-618-2291



info@medtrica.com  
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave  
101-128  
Redondo Beach, CA 90278  
t: 877-703-1079