


Scope Protection




MEDTRICA

TIP PROTECTORS

| Description | Code | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------------------------------------------------|
| <p>Semi-rigid expandable mesh protects the tip/lens of endoscopes from front and side impact. Also for protecting sensitive surgical instruments.</p> <p>Tip Protector also allows moisture to escape, eliminating potential bacteria growth (Pseudomonas species) that can occur with most foam or solid protectors</p> <p>Tip Protector is not abrasive / will not scratch any part of the scope / lens</p> <p>Tip Protector securely fits many different sizes of instruments / does not fall off during transport & storage</p> <p>Options:</p> <p>Small - Red: Protects instruments that are 2mm - 10mm wide - 100 ct / 1 case</p> <p>Large - Blue: Protects instruments that are 7mm - 15mm wide - 100 ct / 1 case</p> | <p>TIP550100</p> <p>TIP650100</p> |  |

FLEXIBLE SCOPE SLEEVE

| | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Paper based scope sleeve protects and wicks moisture away from the insertion tube. It doesn't allow moisture to build up that can result in bacteria growth in-between uses - (Pseudomonas species)</p> <p>Provides surface for identification label placement</p> <p>Options:</p> <p>Flexible Scope Sleeve - Small - 200 - 36" sleeves</p> <p>Flexible Scope Sleeve - Large - 100 - 78" sleeves</p> | <p>ETF500100</p> <p>ETT500050</p> |  <p>Small Scope Sleeve</p> <p>Large Scope Sleeve</p> |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|

Manufactured by:

MEDTRICA SOLUTIONS LTD.
104A-6350 204th Street
Langley, BC Canada V2Y 2V1
t: 604-618-2291



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
101-128
Redondo Beach, CA 90278
t: 877-703-1079