



Medtrica Tip Protector or Foam Protector



<i>Tip Protector Expandable Mesh</i>	<i>Foam Tip Protector</i>
	
Infection Control	
<i>Non-Porous</i>	<i>Porous</i>
<i>Non-Absorbent</i>	<i>Absorbent</i>
<i>Allows Unrestricted Drainage</i>	<i>Restricts Drainage</i>
<i>Allows Unrestricted Air Flow</i>	<i>Restricted Air Flow</i>
Details	
<i>2 Sizes fit 2mm – 15mm</i>	<i>Fits 7mm – 19mm</i>
<i>Secure Fit – will not fall off</i>	<i>Insecure fit – will slip</i>
<i>Protects from Side and Front impact</i>	<i>Protects from Side impact</i>
<i>Use during Storage and all Transport</i>	<i>Use during Transport after a procedure as 'bumper guard'</i>