

Medtrica
4001 Inglewood Ave; 101-128,
Redondo Beach, CA 90278
P: 877-703-1079/F: 310-263-9098
www.medtrica.com



Re the Pressure Mapping test numbers:

What is important to note is that the dividing line between pressure relief and pressure reduction/redistribution is approx 33 mm/Hg.

The most important number to look at is the Average.

Any number **under** 33 mm/Hg shows that the surface is a pressure relief surface for that specific test - (whether it be heel, scapula, etc) - Any number **over** 33 mm/Hg shows that the surface is a pressure reduction surface for that specific test. The higher the number is beyond 33 mm/Hg tells you that the pressure reduction qualities are diminishing.

Please note that most of the results of the testing show that Medtrica's mattresses meet pressure **relief** standards even though they are classified as pressure reduction/redistribution surfaces.

In short, quality pressure reduction mattresses will have Average numbers as close to or below 33mm/Hg

Please also note the length of warranty and the weight ratings on the mattresses, they are well above industry standards.

Best Regards

Dan Doubroff
National Sales Director