

Medtrica Stretcher Mattress Models

	Standard w/ Vinyl	Standard w/ 4 way Stretch	UltraFlex w/ 4 way Stretch
			
Construction	Single Layer	Single Layer	Triple Layer
Memory Foam Layer	No	No	Yes - 4.0 lb Visco elastic
Warranty	2 year Non Prorated on Cover and Core	1 year Non Prorated on Cover and Core	2 year Non Prorated on Cover and Core
Liner	Optional	Optional	Yes
Low Slip Bottom	Yes	Yes	Yes
Adhesive Velcro	Yes	Yes	Yes
Fire Rating	All materials used exceed flammability standards CalTB117, 16CFR1632	All materials used exceed flammability standards CalTB117, 16CFR1632	All materials used exceed flammability standards CalTB117, 16CFR1632
Anti - Bacterial, Anti - Fungal, Anti - Static. Stain & Tear resistant. Latex Free	Yes	Yes	Yes
100% Medical Grade Foam & Fabric	Yes	Yes	Yes
Options	Square or Tapered Corners	Square or Tapered Corners	Square or Tapered Corners
Standard Sizes	Lengths 75" to 78" Widths 24" to 31" Thicknesses 3", 4" & 5"	Lengths 75" to 78" Widths 24" to 31" Thicknesses 3", 4" & 5"	Lengths 75" to 78" Widths 24" to 31" Thicknesses 3", 4" & 5"
	Custom Sizes and Designs Available		
info@medtrica.com		www.medtrica.com	

Medtrica Stretcher Mattress Models

Medtrica offers a full range of 100% medical grade foam stretcher mattresses.

From the single layer Standard Stretcher Mattress to the 3 layer UltraFlex Stretcher Mattress, Medtrica has a stretcher mattress that is the perfect solution for your requirements regarding:

- Patient Safety
- Patient Comfort
- Facility Budget



Medtrica stretcher mattresses provide optimum value because they are engineered for maximum patient safety and comfort while easily fitting into the common budget concerns of your facility.

100% medical grade foam of varying Densities & IFD are used for the construction of each foam layer.

2 fluid proof cover fabric options are available.


4 way Stretchable Pressure Reducing Fabrics are the world leader in Healthcare cover fabrics. They are engineered to provide the best in MVT, pressure management abrasion and tear resistance while maintain optimal performance throughout the entire lifespan of the stretcher mattress.

Anti-Shear liners create less friction and heat buildup between the user and the mattress thus creating a more comfortable, cooler and safer surface.

Anti-Shear liners come standard with UltraFlex models and are optional with Standard models.

All stretcher mattresses meet and exceed OEM specifications.

All Standard and Custom Sizes Available to fit all makes and models of stretcher frames.



Full Non-prorated warranties on Cover and Core Custom Sizing and Construction Available



Manufactured by:

MEDTRICA SOLUTIONS LTD

3525 191 St, Unit 1
Surrey, BC Canada V3Z 0P6

t: 877-703-1079

f: 310-263-9098



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave 101-128
Redondo Beach, CA 90278

t: 877-703-1079

f: 310-263-9098