

# UltraFlex Stretcher Mattress



## Top Layer:

Constructed of 4.0 lb  
Viscoelastic Memory foam  
• 1.5" thick

## Fluid Proof Cover:

4 way Stretch Pressure Reducing

## Anti-Shear Liner:

Dissipates Heat & Reduces Friction

## Middle Layer:

Constructed of  
1545 lb density foam  
• 1.5" or 2.5" thick

## Bottom Layer:

Constructed of  
1890 lb density foam  
• 1" thick

**AVAILABLE SIZES:**  
All standard & Custom  
sizes Available to Fit all  
Makes and Models of  
Stretcher Frames

## GUARANTEE:

A full non-prorated  
guarantee on both the  
cover and the core for all  
stretcher mattresses.

Adhesive Velcro  
included to keep the  
mattress secure to  
the stretcher.

## FLAMMABILITY

All materials meet Cal  
TB117 and CFR1632.  
Option to meet CalTB129  
& Boston FC IX 11.

All Stretcher mattresses  
meet & exceed OEM  
specifications.

## CHOICE OF MATTRESS CONSTRUCTION:

- 4" or 5" Thick; All standard and custom Lengths & Widths
- Tapered or Square corners.

## COVER:

The Pressure Reducing stretchable cover fabric is 100% medical grade, anti-fungal, anti-bacterial, latex-free and completely fluid-proof. The UltraFlex model comes standard with an inner Anti-Shear Liner which dissipates heat and reduces friction for greater comfort and safety

## 4 WAY STRETCHABLE PRESSURE REDUCING FABRIC:

4 way Stretchable Pressure Reducing Fabric is the world leader in Healthcare cover fabric. The fabric is engineered to provide the best in MVT, pressure management, abrasion & tear resistance while maintaining optimal performance throughout the entire lifespan of the stretcher mattress. Latex Free.

## FOAM CORE:

3 layers of 100% medical grade non-recycled, non-allergenic, HD polyurethane & Viscoelastic (memory) foam without fillers or additives. Latex Free.

## CLEANING:

Stretcher mattress is disinfected and wiped clean using standard cleaning products.

## MORE DURABLE & VERSATILE

*Medtrica stretcher mattresses  
provide greater comfort &  
durability because of  
unsurpassed Quality.*

*100% medical grade foam of  
varying Densities & IFD are used  
for the construction of each  
foam layer.*

*Manufactured in the  
USA & Canada -  
Engineered to Last.*

[info@medtrica.com](mailto:info@medtrica.com)  
[www.medtrica.com](http://www.medtrica.com)



*A better quality stretcher mattress, more durable - for a better price.*