




PROTECTION

EndoCap™ Tip Protector


- Semi-rigid expandable mesh protects the tip/lens of endoscopes from front and side impact.
- Also for sensitive surgical instruments.
- Allows moisture to escape, eliminating potential bacteria growth
- Will not fall off and is not abrasive/will not scratch any part of the scope/lens
- For use during transport and storage
- Securely fits many different sizes of instruments
- Each case contains 10 packs of 10 (clean) or 50 units (sterile)

Clean or Sterile

Part #	Short Description	Quantity / Pack Size	Product Description	Image
EndoCap™ Small Tip Protector				
TIP550100	TIP / LENS PROTECTOR – SMALL	100 ct/case (10 packs of 10)	Tip Protection for Scope Lens - Small Red • Protects instruments that are 2mm - 10mm wide	
EndoCap™ Large Tip Protector				
TIP650100	TIP / LENS PROTECTOR – LARGE	100 ct/case (10 packs of 10)	Tip Protector for Scope Lens - Large Blue • Protects instruments that are 7mm - 15mm wide	
EndoCap™ STERILE Small Tip Protector				
TIP550050S	STERILE TIP / LENS PROTECTOR – SMALL	50 ct/case 50 packs of 1 ea	Sterile Tip Protection for Scope Lens - Small Red • Protects instruments that are 2mm - 10mm wide	
EndoCap™ STERILE Large Tip Protector				
TIP650050S	STERILE TIP / LENS PROTECTOR – LARGE	50 ct/case 50 packs of 1 ea	Sterile Tip Protector for Scope Lens - Large Blue • Protects instruments that are 7mm - 15mm wide	

EndoSock™ Scope Sleeve

- Provides a barrier from Physical contact
- Provides a barrier filter from the Environment
- Does not allow moisture build up and the resulting potential bacteria growth from dampness on the insertion tube / wicks moisture away
- Provides a surface for the placement of labels showing detailed scope cleaning information for reference, without having to place the label on the scope body (labels can be easily removed from the EndoSock and placed on a patients chart, etc)
- Does not need to be taped directly onto the endoscope boot

Part #	Short Description	Quantity / Pack Size	Product Description	Image
Endosock Small Scope Sleeve				
EPS369502	SCOPE SLEEVE - SMALL - 36"	200 ct/case (10 packs of 20)	Flexible Scope Sleeve: Small - 36" Long - 200 count	
Endosock Large Scope Sleeve				
EPS789501	SCOPE SLEEVE - LARGE - 78"	100 ct/case (10 packs of 10)	Flexible Scope Sleeve: Large - 78" Long - 100 count	